

Medicinal Stain Remover

Medicinal Stain Remover is highly effective for removing a wide range of stains caused by medicinal products. **Medicinal Stain Remover** is a blend of active solvents and surfactants. It is able to remove difficult stains without harming cotton and synthetic fibers. Phosphate Free. For institutional laundry use only.

- Solvent blend
- Contains surfactants
- Ready-to-use formula

Technical Specifications

Appearance	Clear liquid
Dilution Rate	1:1; 4-5 OZ/GAL
Color	Colorless
Scent	Solvent
Non Volatile Matter	12.5 - 13.3%
pH @ 25 deg. C	1% 11.3 - 11.9



Safety

See material safety data sheet and product for safety information, handling and proper use.

Availability

HIL0005706

4 - 1 Gallon Containers

Directions

PRE-WASH STAIN TREATMENT:

Dilute **MEDICINAL STAIN REMOVER** with water 1:1. Spray directly on fabric stain. Wash fabric through normal cycle.

POST-WASH STAIN TREATMENT:

Add 4-5 ounces of **MEDICINAL STAIN REMOVER** per gallon of hot water (140-160°F). Allow stained fabric to soak in solution for 20-30 minutes. Wash fabric through normal cycles, do not rinse before laundering.

PRECAUTION:

For colored fabric or unknown fiber content, test on hidden area prior to use.